

E-Tank Duration Chart

Guide for Estimating Minutes of Available Oxygen

Liters per Minute	Pressure (PSI) in E cylinder							
	500	750	1000	1250	1500	1750	2000	2200
	Approximate Minutes of Operation							
1	140	210	280	350	420	490	560	616
2	70	105	140	175	210	245	280	308
3	46	70	93	116	140	163	186	205
4	35	52	70	87	105	122	140	154
5	28	42	56	70	84	98	112	123
6	23	35	46	58	70	81	93	102
10	14	21	28	35	42	49	56	61
12	11	17	23	29	35	40	46	51
15	09	14	18	23	28	32	37	41

Note: Shaded area indicates “at-risk” minutes. When psi regulator readings are between the psi readings on the table, use the lower psi reading to determine the approximate minutes of oxygen available.



**Maverick
Oxygen &
Respiratory
Equipment, LLC**